



Sample Nanomanipulation of 1-D Nanostructures

Using attocube's ECSx3030 Positioners Inside a Scanning Electron Microscope

Rodrigo Bernal and Horacio Espinos
Mechanical Engineering Department, Northwestern University,
Evanston IL, USA

Introduction

The characterization of 1-D nanostructures such as nanowires and nanotubes has recently become a topic of major interest, due to the potential of these nanostructures to be employed in the next generation of advanced materials and electronic devices. Due to their small characteristic size (< 100 nm in diameter) and an increased surface to volume ratio, many nanostructures of a given material display significantly enhanced properties compared to the macroscale, bulk material. For example, zinc oxide (ZnO) nanowires, display a greater modulus of elasticity than bulk for diameters less than ~ 80 nm [1]. Similar behavior has also been discovered for other semiconducting nanowires such as gallium nitride (GaN) [2], and metallic nanowires such as silver [3].




Figure 1: MEMS device for the characterization of electromechanical properties of nanowires.

The emergence of these size-effects in the properties of 1-D nanostructures has therefore elicited a need for the characterization and unambiguous measurement of these properties, as they are critical for the development, design and robustness of future applications employing nanostructures as their functional elements. However, the small size of the specimens imposes significant challenges for specimen preparation and testing. To overcome these difficulties, Prof. Horacio Espinosa's group at the Mechanical Engineering Department in Northwestern University, USA, has developed a microelectromechanical system (MEMS) that allows uniaxial mechanical testing of nanowires, with nanometer (nm) and nanoNewton (nN) resolution, therefore allowing the accurate measurement of mechanical and failure properties such as the elastic modulus, yield and fracture strengths [4] (see Figure 1). The system can also be employed to carry out simultaneous four-point electrical measurements as the mechanical straining is taking place, therefore allowing measurements of

electromechanical properties of nanowires, such as piezo-resistivity and piezo-electricity [5].

Setup

A critical step in these experiments is the sample nanomanipulation, where a nanowire must be transferred to the testing MEMS. This delicate procedure requires an instrument with capabilities of nm resolution in displacement, smooth movement to avoid any vibrations that may harm the sample and seamless integration with an scanning electron microscope (SEM), to allow observation of the preparation process (an optical microscope lacks the needed resolution). A large range of total movement (> 1 mm) is also convenient to allow flexibility in the setup of the experiment. In recent experiments, the Espinosa group has employed an attocube nanomanipulator, composed of three stacked ECSx3030 positioners, one for each axis of movement, in order to accomplish this task. The nanomanipulator is positioned inside the SEM chamber and interfaced to one of attocube's piezo motion controllers, located outside the chamber, through vacuum feedthroughs. The controller is connected via USB to a laptop furnished with attocube software for controlling the manipulator.




Figure 2: Silver nanowire (yellow) laying on a copper TEM grid (blue) ready to be manipulated. The tungsten tip (red) attached to the nanomanipulator is used to approach finely to the specimen.

Measurement Results

To perform the nanomanipulation, a sharp tungsten tip is attached to the manipulator (Figure 2, manually colored). Coarse-steps, of the order of μm , are used to approach the nanowire with the tip. Afterwards, fine steps with magnitudes



Sample Nanomanipulation of 1-D Nanostructures

Using attocube's ECSx3030 Positioners Inside a Scanning Electron Microscope




Figure 3: Silver nanowire (yellow) picked from the grid and ready to be mounted on the MEMS device of Figure 1.

microscopy testing of nanostructures. *Journal of Microelectromechanical Systems* 16, 1219 (2007).

- [5] R.A. Bernal, T. Filleter, J.G. Connell, K. Sohn, J. Huang, L.J. Lauhon, and H.D. Espinosa, In Situ Electron Microscopy Four-Point Electromechanical Characterization of Freestanding Metallic and Semiconducting Nanowires. *Small* 10, 725 (2014).

H.E. acknowledges funding from the US Army Research Office through DURIP award No. W911NF-12-1-0366.

controllable from sub-100 nm to few nm, are used to establish gentle contact between the tip and nanowire.

Once the nanowire is picked (Figure 3, manually colored), a similar approach is used to position it on top of the MEMS device. Final attachment of the specimen is carried out by electron beam induced deposition of platinum (EBID-Pt) after which the tip is retracted from the area of interest.

Conclusion

This Application Note shows the applicability of attocube's positioners to carry out nanometer precise movements without inducing unwanted interactions in a complex setup such as an SEM. The ECSx3030 is perfectly suitable as it can be stacked easily to a multi-axis solution with the help of some adapter plates and/or L-holder, while still matching the requirements for precision.

References

- [1] R. Agrawal, B. Peng, E.E. Gdoutos, and H.D. Espinosa, Elasticity size effects in ZnO nanowires - A combined Experimental-Computational approach. *Nano Letters* 8, 3668 (2008).
- [2] R.A. Bernal, R. Agrawal, B. Peng, K.A. Bertness, N.A. Sanford, A.V. Davydov, and H.D. Espinosa, Effect of Growth Orientation and Diameter on the Elasticity of GaN Nanowires. A Combined in Situ TEM and Atomistic Modeling Investigation. *Nano Letters* 11, 548 (2011).
- [3] T. Filleter, S. Ryu, K. Kang, J. Yin, R. Bernal, K. Sohn, S. Li, J. Huang, W. Cai, and H.D. Espinosa, Nucleation-Controlled Distributed Plasticity in Penta-Twinned Silver Nanowires. *Small* 8, 2986 (2012).
- [4] H.D. Espinosa, Y. Zhu, and N. Moldovan, Design and operation of a MEMS-based material testing system for in-situ electron