

LT-APO/ULWD/VISIR/0.35/xs

1014253

Technical Specifications

Optical Data	
clear aperture [mm]	4.7
focal length [mm]	6.71
numerical aperture(NA)	0.35
working distance	12.0 mm
Spectral Performance	
AR coating (> 80% transmission) [nm]	400 .. 1000
apochromatic range (df < +/- delta) [nm]	520 .. 870
Compatibility	
environment	high magnetic fields
	low temperature
	ultra-high vacuum
compatible setups	attoCFM I/cust, attoDRY800, AFM/CFM/cust
Size and Dimensions	
length	50.35 mm
weight	approx. 50 g



δf : chromatic focal shift, $\Delta = n \cdot \lambda_{\text{ref}} / (2 \cdot \text{NA}^2)$: depth of focus, n : refractive index, λ_{ref} : wavelength used to define focal plane with max. Δ



