

# LT-LWDO/0.55 low temperature compatible aspheric objective for attoCFM

I,  
1007212

## Technical Specifications

Optical Data	
clear aperture [mm]	4.95
focal length [mm]	4.51
numerical aperture(NA)	0.55
working distance	2.7 mm
Spectral Performance	
AR coating (> 80% trasnmission) [nm]	A : 350 .. 700
	B : 600 .. 1050
	C : 1050 .. 1700
Compatibility	
environment	high magnetic fields
	high vacuum
	low temperature
compatible setups	attoCFM I, attoCFM IV, attoDRY800
mounting	M8 x 0.5
	metric M12 x 0.5
Size and Dimensions	
length	31.8 mm



