

LT-IWDO/0.68 low temperature compatible aspheric objective for attoCFM

I,
1009849

Technical Specifications

Optical Data	
clear aperture [mm]	5
focal length [mm]	3.1
numerical aperture(NA)	0.68
working distance	1.5 mm
Spectral Performance	
AR coating (> 80% trasnmission) [nm]	A : 350 .. 700
	B : 600 .. 1050
	C : 1050 .. 1700
Compatibility	
environment	high magnetic fields
	high vacuum
	low temperature
compatible setups	attoCFM I, attoCFM IV, attoDRY800
mounting	M8 x 0.5
	metric M12 x 0.5
Size and Dimensions	
length	32.9 mm



