

# AMCC/1/RES/5-pin/2.0/RT - connection cable

1020988

## Technical Specifications

General Specifications	
connector (controller side "1")	D-Sub HD 26-pin (male)
connector (feedthrough or positioner side "2")	5-pin socket (female)
bending diameter (static)	= OD * 5 [mm]
bending diameter (dynamic)	= OD * 10 [mm]
Size and Dimensions	
cable length	2 m
cable diameter OD	8.7 mm
weight	152 g
Working Conditions	
operating pressure	ambient
operating temperature	ambient
Materials	
body connector (controller side "1")	die-cast zinc Alloy, tin, copper alloy; gold-plated, glass filled thermoplasti
body connector (feedthrough or positioner side "2")	PA 4.6. GF, CuZn-alloy, tin, gold, CuBe-alloy, Ni+0,25µm Au, Ni+4...6µm Sn
cable wire	copper
cable isolation	PVC (outer), PE (inner)
Scope of Delivery	
cable	AMCC/RES cable to connect controller with AN* Positioner

