

# AMCC/1/NUM/14-Binder/2.0/RT - connection cable

1020999

## Technical Specifications

General Specifications	
connector (controller side "1")	D-Sub HD 26-pin (male)
connector (feedthrough or positioner side "2")	14 pin circular binder (female)
bending diameter (static)	= OD * 5 [mm]
bending diameter (dynamic)	= OD * 10 [mm]
Size and Dimensions	
cable length	2 m
cable diameter OD	7.9 mm
weight	215 g
Working Conditions	
operating pressure	ambient
operating temperature	ambient
Materials	
body connector (controller side "1")	die-cast zinc Alloy, tin, copper alloy; gold-plated, glass filled thermoplasti
body connector (feedthrough or positioner side "2")	PA 66 (UL 94 HB), copper tin, gold-plated, PBT (UL 94 V-0)
cable wire	copper
cable isolation	PVC (outer), PE (inner)
Scope of Delivery	
cable	AMCC/NUM cable to connect controller with EVFT

