














Accessories

Sample Holders

upgrade options for microscopes and/or cryostats

	Sample holder	Base & socket materials	quick exchange	# of electrical contacts	30 cm wires with connector incl.	type of wires	heater & calibrated temp. Sensor	He exchange gas (attoDRY1000/1100/2100)	vacuum (attoDRY 800)	footprint	spare sample holder	Art. No.
	ASH/QE/0	Ti base, Ti slider	•	0	n.a.	n.a.	•	•		101/102	1006288	1006287
	ASH/QE/4CX	Ti base, PCB holder	•	8	•	4 coaxials	•	•		101/102	1009161 (thick) 1009187 (thin)	1013376
	Tray for 3x ASH/QE/4CX sample holders		•	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.		1011525
	ASH/QE/8/CFM	Ti base, PCB holder	•	8	•	8 wires as 4 twisted pairs	•	•		101/102	1012803	1012172
	ASH/PCB/12	Ti base, PCB holder		12	•	12 wires as 6 twisted pairs		•	•	101/102	1005708	1005710
	ASH51/PCB/GND	PCB holder		GND	•	1 coaxial	n.a.	•	•	50/51	1007920	1013377
	ASH51/PCB/4+GND	PCB holder		4+GND	•	4 coaxials	n.a.	•	•	50/51	1008130	1013378
	ASH/SWD			0	n.a.	n.a.	•		•	101/102		1012785
	ASH/Zvar			0	n.a.	n.a.	•		•	attoDRY800 cold plate		1013136
	ASH/PCB/12mini			12	•	12 wires as 6 twisted pairs	n.a.		•	n.a.		1013138
	ASH/CC/20	LCCC	•	20	•	20 wires as 10 twisted pairs	n.a.	•		ANR51	1007718	1007229
	Base plate	Ti base	n.a.	n.a.	n.a.	n.a.	•	•	•	101/102	n.a.	1006900
	Base plate	Cu base	n.a.	n.a.	n.a.	n.a.	•	•	•	50/51	n.a.	1004347