


Accessories

Adapter Plates EAP

EAP adapter plates are used for mounting attocube's ECS positioners with vertical orientation and cross-mounting of two differently sized ECS positioners.



Article	Material	
	/Al	/StSt
EAP1/2 – L-bracket for ECS3030	1005784	1006818
EAP3/4 – L-bracket for ECS5050	1007996	1008002
EAP5/6 – Combination of differently sized ECS	1008959	1008971
EAP7/8 – L-bracket for ECS3040 & ECS3050	1010272	1010274
EAP9/10 – L-bracket for ECS3060, ECS3070 & ECS3080	1010275	1010276
EAP11 – ECR4040 Cage Mount	1010385	
EAP12 – adapter plate for ECR4040	1010386	
Magnetic sample holder	1007589	

Combining Industrial Line Positioners - Overview

ECSx3030	direct	direct	direct	direct	direct	direct	direct	direct	direct	direct	-	-	direct	direct
ECSx3040	direct	direct	direct	direct	direct	direct	direct	direct	direct	direct	-	-	direct	direct
ECSx3050	direct	direct	direct	direct	direct	direct	direct	direct	direct	direct	-	-	direct	direct
ECSx3060	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	-	EAP5/6	EAP5/6	EAP5/6
ECSx3070	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	-	EAP5/6	EAP5/6	EAP5/6
ECSx3080	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	EAP5/6	-	EAP5/6	EAP5/6	EAP5/6
ECSx5050	-	-	-	-	-	-	direct	direct	-	direct	-	-	*	direct
ECSxy5050	-	-	-	-	-	-	direct	direct	-	direct	-	-	*	direct
ECSz3030	direct	direct	direct	direct	direct	direct	direct	direct	direct	direct	-	-	direct	direct
ECSz5050	-	-	-	-	-	-	direct	direct	-	direct	-	-	*	direct
ECR3030	direct	direct	direct	direct	direct	direct	direct	direct	direct	direct	-	-	direct	direct
ECR4040	-	-	-	-	-	-	direct	direct	-	direct	-	-	*	direct
ECR5050	-	-	-	-	-	-	direct	direct	-	direct	-	-	-	direct
ECG5050	-	-	-	-	-	-	direct	direct	-	direct	-	-	*	direct
mount on	ECSx3030	ECSx3040	ECSx3050	ECSx3060	ECSx3070	ECSx3080	ECSx5050	ECSxy5050	ECSz3030	ECSz5050	ECR3030	ECR4040	ECR5050	ECG5050

* on request