

1007211

1007211

Technical Specifications

Optical Data	
clear aperture [mm]	5
focal length [mm]	3.1
numerical aperture(NA)	0.68
working distance	1.5 mm
Spectral Performance	
AR coating (> 80% transmission) [nm]	A : 350 .. 700 B : 600 .. 1050 C : 1050 .. 1700
Compatibility	
environment	high magnetic fields high vacuum low temperature
compatible setups	attoCFM II, attoCFM IV, attoDRY800
mounting	clamped
Size and Dimensions	
length	32.9 mm



