

# 1007212

1007212

## Technical Specifications

Optical Data	
clear aperture [mm]	4.95
focal length [mm]	4.51
numerical aperture(NA)	0.55
working distance	2.7 mm
Spectral Performance	
AR coating (> 80% transmission) [nm]	A : 350 .. 700 B : 600 .. 1050 C : 1050 .. 1700
Compatibility	
environment	high magnetic fields high vacuum low temperature
compatible setups	attoCFM I, attoCFM IV, attoDRY800
suitable broadband collimator	RT-APO/VIS-NIR/0.13
mounting	M8 x 0.5 metric M12 x 0.5
Size and Dimensions	
length	31.8 mm



