

# AMCC/3/RES/D-Sub-15/2.0/(U)HV - connection cable

1020369

## Technical Specifications

General Specifications	
connector (controller side "1")	3 x D-Sub HD 26-pin (male)
connector (feedthrough or positioner side "2")	D-Sub 15-pin (male)
bending diameter (static)	= OD * 5 [mm]
bending diameter (dynamic)	= OD * 10 [mm]
Size and Dimensions	
cable diameter OD	8.7 mm
weight	530 g
Working Conditions	
operating pressure	ambient
operating temperature	ambient
Materials	
body connector (controller side "1")	die-cast zinc Alloy, tin, copper alloy; gold-plated, glass filled thermoplast
body connector (feedthrough or positioner side "2")	die-cast zinc alloy, steel/tin-plated, polybutylenterephthalat PBT, brass: gold-flashed
cable wire	copper
cable isolation	Fluorthermoplast
Scope of Delivery	
cable	AMCC/RES cable to connect 3 controller axis with one EVFT
adapters	SubD15 male-male (feedthrough mockup)

