

ANGp101/LT/HV - goniometer for "phi" tilting

1002742

Technical Specifications

| | |
|--------------------------------------|--------------------------------|
| Size and Dimensions | |
| footprint; height | 24 mm x 24 mm; 11 mm |
| distance center of rotation to top | 51 mm |
| maximum installation space | 28.6 mm x 24 mm; 11.8 mm |
| weight | 18 g |
| Materials | |
| positioner body | titanium |
| actuator | PZT ceramics |
| connecting wires | insulated twisted pair, copper |
| Coarse Positioning Mode | |
| travel range (step mode) | 5.4 ° |
| maximum drive velocity @ 300 K | ~ 1 °/s |
| input voltage range | 0 - 60 V |
| Fine Positioning Mode | |
| fine positioning resolution | no fine positioning capability |
| Load (@ ambient conditions) | |
| maximum load | 100 g |
| maximum dynamic force along the axis | 2 N |

