

# ANGt101/LT/UHV - goniometer for "theta" tilting

1002737

## Technical Specifications

<b>Size and Dimensions</b>	
footprint; height	24 mm x 24 mm; 11 mm
distance center of rotation to top	40 mm
maximum installation space	28.6 mm x 24 mm; 11.8 mm
weight	18 g
height	11 mm
<b>Materials</b>	
positioner body	titanium
actuator	PZT ceramics
connecting wires	insulated twisted pair, copper
<b>Coarse Positioning Mode</b>	
travel range (step mode)	6.6 °
maximum drive velocity @ 300 K	~ 1 °/s
input voltage range	0 - 60 V
<b>Fine Positioning Mode</b>	
fine positioning resolution	no fine positioning capability
<b>Load (@ ambient conditions)</b>	
maximum load	100 g
maximum dynamic force along the axis	2 N
<b>General Specifications</b>	
environment	/LT/UHV

