

# ANPx311/RES/HV - linear x-nanopositioner

1008591

## Technical Specifications

<b>Size and Dimensions</b>	
footprint; height	34 mm x 30 mm; 10 mm
maximum installation space	34 mm x 36 mm; 10 mm
weight	32 g
<b>Materials</b>	
positioner body	titanium
actuator	PZT ceramics
connecting wires	insulated twisted pair, copper
<b>Coarse Positioning Mode</b>	
travel range (step mode)	6 mm
maximum drive velocity @ 300 K	~ 3 mm/s
<b>Fine Positioning Mode</b>	
fine linear positioning range @ 300 K	5 $\mu$ m
fine positioning resolution	sub-nm
input DC voltage range @ 300 K	0 - 100 V
<b>Position Encoder</b>	
readout mechanism	resistive sensor
sensor resolution	~ 200 nm
repeatability	1..2 $\mu$ m (unidirectional)
<b>Load (@ ambient conditions)</b>	
maximum load	2 kg
maximum dynamic force along the axis	2 N
<b>General Specifications</b>	
environment	/HV
min. first natural frequency	2600 Hz

