

ANPz51/HV - linear z-nanopositioner

1002225

Technical Specifications

Size and Dimensions	
footprint; height	15 mm x 15 mm; 13.5 mm
maximum installation space	15 mm x 15 mm; 16 mm
weight	11 g
Materials	
positioner body	titanium
actuator	PZT ceramics
connecting wires	insulated twisted pair, copper
Coarse Positioning Mode	
travel range (step mode)	2.5 mm
maximum drive velocity @ 300 K	~ 1 mm/s
Fine Positioning Mode	
fine linear positioning range @ 300 K	3.5 μ m
fine positioning resolution	sub-nm
Load (@ ambient conditions)	
maximum load	50 g
maximum dynamic force along the axis	1 N

