

# ANR101/UHV - rotator (360° endless)

1002205

## Technical Specifications

<b>Size and Dimensions</b>	
footprint; height	24 mm x 24 mm; 15.2 mm
weight	36 g
height	15.2 mm
<b>Materials</b>	
positioner body	titanium
actuator	PZT ceramics
connecting wires	insulated twisted pair, copper
<b>Coarse Positioning Mode</b>	
travel range (step mode)	360 ° endless
maximum drive velocity @ 300 K	~ 30 °/s
input voltage range	0 - 60 V
typical actuator capacitance @ 300 K [µF]	1.2
<b>Fine Positioning Mode</b>	
fine angular positioning range @ 300 K	70 m°
fine positioning resolution	µ°
input DC voltage range @ 300 K	0 - 120 V
<b>Load (@ ambient conditions)</b>	
maximum load	100 g
maximum dynamic torque around axis	0.8 Ncm
<b>General Specifications</b>	
environment	/UHV

