

# ANRv220/RES/LT - rotator (360° endless)

1006366

## Technical Specifications

<b>Size and Dimensions</b>	
footprint; height	27 mm x 12 mm; 27 mm
weight	38 g
<b>Materials</b>	
positioner body	beryllium copper
actuator	PZT ceramics
connecting wires	insulated twisted pair, copper
<b>Coarse Positioning Mode</b>	
travel range (step mode)	360 ° endless
maximum drive velocity @ 300 K	~ 30 °/s
<b>Fine Positioning Mode</b>	
fine angular positioning range @ 300 K	70 m°
fine angular positioning range @ 4 K	14 m°
fine positioning resolution	μ°
<b>Position Encoder</b>	
readout mechanism	resistive sensor
encoded travel range	315°
sensor resolution	~ 6 m°
repeatability	50 m° (unidirectional)
<b>Load (@ ambient conditions)</b>	
maximum load	100 g
maximum dynamic torque around axis	1 Ncm

