

# ANSx150/LT/UHV x-Scanner made from Titanium

1006873

## Technical Specifications

<b>Size and Dimensions</b>	
footprint; height	24 mm x 24 mm; 9 mm
maximum installation space	24 mm x 24 mm; 9 mm
weight	18 g
<b>Materials</b>	
positioner body	titanium
actuator	PZT ceramics
connecting wires	insulated twisted pair, copper
<b>Coarse Positioning Mode</b>	
input voltage range	0 - 60 V
<b>Fine Positioning Mode</b>	
fine linear positioning range @ 300 K	80 $\mu$ m
fine linear positioning range @ 4 K	125 $\mu$ m
fine positioning resolution	sub-nm
typical actuator capacitance fine @ 300 K	11.1 $\mu$ F
typical actuator capacitance fine @ 4 K	2.6 $\mu$ F
input DC voltage range @ 300 K	0 - 60 V
input DC voltage range @ 4 K	+/- 120V
<b>Load (@ ambient conditions)</b>	
maximum load	100 g

