

ANSxy50/LT/UHV xy-Scanner made from Titanium

1001089

Technical Specifications

Size and Dimensions	
footprint; height	15 mm x 15 mm; 7 mm
maximum installation space	15 mm x 15 mm; 7 mm
weight	7 g
Materials	
positioner body	titanium
actuator	PZT ceramics
connecting wires	insulated twisted pair, copper
Fine Positioning Mode	
fine linear positioning range @ 300 K	30 x 30 μm
fine linear positioning range @ 4 K	15 x 15 μm
fine positioning resolution	sub-nm
typical actuator capacitance fine @ 300 K	1.1 μF
typical actuator capacitance fine @ 4 K	0.2 μF
input DC voltage range @ 300 K	0 - 60 V
input DC voltage range @ 4 K	0 - 150 V
*	scan range depends on environment
Load (@ ambient conditions)	
maximum load	50 g
General Specifications	
environment	/LT/UHV

