

IGLU Compressor

1023623

Technical Specifications

System Description	
helium compressor for closed-cycle cryocoolers	.
air cooling	speed regulation - „adaptive power“
orientation free operation	ethernet interface for remote control
General Specifications	
mounting	19 inch rack mounting with removeable wheels and handles
pressure range	0..25 bar / 0..362 psi
ambient temperature	10..30 °C (non-condensing)
dimensions	446 x 266 x 641 mm, 17.5" x 10.5" x 25.25"; excluding handles, wheels and front panel
weight	approx. 70 kg
country of origin	Federal Republic of Germany
Electrical Specifications	
maximum input power	< 1.4 kW
AC input voltage and line frequency	single phase 230 V, 50 Hz/60 Hz 110 V option with external transformer
inverter driven pump incl. EMC filter	.

