

LT-APO/LWD/NIR/0.63: low temperature compatible, apochromatic objective

1013375

Technical Specifications

Optical Data	
clear aperture [mm]	4.7
focal length [mm]	3.73
numerical aperture(NA)	0.63
working distance	5.0 mm
Spectral Performance	
apochromatic range (df < +/- delta) [nm]	725 .. 970
Compatibility	
environment	high magnetic fields high vacuum low temperature
compatible setups	attoCFM I/cust, attoDRY800, AFM/CFM/cust
Size and Dimensions	
length	81 mm
weight	approx. 100 g



δf : chromatic focal shift, $\Delta = n \cdot \lambda_{ref} / (2 \cdot NA^2)$: depth of focus, n : refractive index, λ_{ref} : wavelength used to define focal plane with max. Δ



