

LT-APO/VIS/0.82/xs

1011768

Technical Specifications

Optical Data	
clear aperture [mm]	4.7
focal length [mm]	2.87
numerical aperture (NA)	0.82
working distance	0.65 mm (1.40 mm)
Spectral Performance	
AR coating (> 80% transmission) [nm]	400 .. 1000
apochromatic range (df < +/- delta) [nm]	465 .. 590
Compatibility	
environment	high magnetic fields high vacuum low temperature
compatible setups	attoCFM Ixs
suitable broadband collimator	RT-APO/VIS-NIR/0.13
Size and Dimensions	
length	50.35 mm
weight	43 g



δf : chromatic focal shift, $\Delta = n \cdot \lambda_{ref} / (2 \cdot NA^2)$: depth of focus, n : refractive index, λ_{ref} : wavelength used to define focal plane with max. Δ



