

# LT-IWDO/0.68 low temperature compatible fiber-based aspheric objective, 1007211

## Technical Specifications

Optical Data	
clear aperture [mm]	5
focal length [mm]	3.1
numerical aperture(NA)	0.68
working distance	1.5 mm
Spectral Performance	
AR coating (> 80% transmission) [nm]	C : 1050 .. 1700
	A : 350 .. 700
	B : 600 .. 1050
Compatibility	
environment	high vacuum
	low temperature
	high magnetic fields
compatible setups	attoCFM II, attoCFM IV, attoDRY800
mounting	clamped
Size and Dimensions	
diameter	12 mm
length	36.4 mm



