

LT-IWDO/0.68 low temperature compatible fiber-based aspheric objective, 1007211

Technical Specifications

| | |
|----------------------------|--|
| Optical Data | |
| clear aperture [mm] | 5 |
| focal length [mm] | 3.1 |
| numerical aperture(NA) | 0.68 |
| working distance | 1.5 mm |
| Compatibility | |
| environment | high magnetic fields high vacuum low temperature |
| compatible setups | attoCFM II, attoCFM IV, attoDRY800 |
| mounting | clamped |
| Size and Dimensions | |
| length | 36.4 mm |



