

Sensor Head D12/F2.8

1011404

Technical Specifications

Sensor Head	
optics type	focusing
dimensions	Ø 12 mm, length 32.3 mm
weight [g]	8.8
mounting	clamped
spot size	14 µm
focal length [mm]	2.8
typical beam diameter at focal point [µm]	~ 14
connector	FC/PC
working range (target: glass) [mm]	2.8 (± 4 µm)
angular tolerance (target: glass) [°]	± 4
working range (target: wafer) [mm]	2.8 (± 14 µm)
angular tolerance (target: wafer) [°]	± 7
working range (target: mirror) [mm]	2.8 (± 20 µm)
angular tolerance (target: mirror) [°]	± 10
fiber connection	screwed
Environment	
pressure	ambient pressure
temperature	0 .. 50 °C
Main Materials	
housing	Titanium Grade 2 (3.7035)
lens	glass (with AR C coating)
Scope of Delivery	
fiber IDS-sensor head	IDSMF/PC/RT (1009980)
number of axes	1
mirror	Ø 9 mm, Al coated (1010846)



