

Toolboxes

1022272

Technical Specifications

Nanopositioner Toolbox "PEEK" (1004256, 1008438)	
Slotted cheese head screws made from titan	M1.6x2 [4 pcs.]
	M1.6x4 [12 pcs.]
	M1.6x6 [4 pcs.]
	M2.5x10 [8 pcs.]
	M2x10 [4 pcs.]
	M2x3 [4 pcs.]
	M2x5 [12 pcs.]
	M2x7 [4 pcs.]
	M2x8 [8 pcs.]
PEEK connectors	2 pin [3 pcs.]
	3 pin [3 pcs.]
	5 pin [3 pcs.]
Tools	Screwdriver blade 1.5 mm
	Screwdriver blade 2.3 mm
	Tweezers
Wire	Copper wire with polyimid coating (32 m x 0.2 mm) [2 pcs.]
Article numbers (adapter plates)	1004256: with adapter plate ANB100 (for AN*5*,10*,220)
	1008438: with adapter plate ANB200 (for all AN* except ANPx(y)312)
Nanopositioner Toolbox "mK" (1022161)	
Slotted cheese head screws made from brass	M1.6x4 [12 pcs.]
	M1.6x6 [4 pcs.]
	M2.5x10 [4 pcs.]
	M2x10 [4 pcs.]
	M2x3 [4 pcs.]
	M2x5 [12 pcs.]
	M2x8 [4 pcs.]
PEEK connectors	2 pin [3 pcs.]
	3 pin [3 pcs.]
	5 pin [3 pcs.]
Tools	Screwdriver blade 1.5 mm
	Screwdriver blade 2.3 mm
	Tweezers
Wire	Copper wire with polyimid coating (32 m x 0.2 mm) [2 pcs.]
Toolbox RES cabling for AMC300 (1021013)	
Cable	1x(2xAWG26) + 1x(2xAWG28) + 5x(2xAWG30) [20 m]
Connector	D-Sub, size A, 26 pin connector (housing and plug) [3 sets]
Equipment	Insulating caps (3 pcs.)
Nanopositioner Toolbox "PA" (1004124, 1008437)	
Slotted cheese head screws made from titan	M1.6x2 [4 pcs.]
	M1.6x4 [12 pcs.]
	M1.6x6 [4 pcs.]
	M2.5x10 [8 pcs.]
	M2x10 [4 pcs.]
	M2x3 [4 pcs.]
	M2x5 [12 pcs.]
	M2x7 [4 pcs.]
	M2x8 [8 pcs.]
PA connectors	10 pin [10 pcs.]
	2 pin [10 pcs.]
	3 pin [10 pcs.]
	5 pin [10 pcs.]
	6 pin [10 pcs.]
Tools	Screwdriver blade 1.5 mm
	Screwdriver blade 2.3 mm
	Tweezers
Wire	Copper wire with polyimid coating (32 m x 0.2 mm) [2 pcs.]
Article numbers (adapter plates)	1004124: with adapter plate ANB100 (for AN*5*,10*,220)
	1008437: with adapter plate ANB200 (for all AN* except ANPx(y)312)

